


# Emergency Disconnect Procedure

Model AK 98  
For Patients  
Graft or Fistula

**1**




**PRESS BLOOD PUMP BUTTON.**  
Blood pump will stop and button will start flashing.

**2**




**USING THE PRODUCTS FOUND IN YOUR EMERGENCY KIT...**

**3**



**CLOSE TWO CLAMPS ON NEEDLES AND BLOODLINES. ADD SAFETY CLAMPS BELOW AND ABOVE THE CONNECTION.**

**4**



**DISCONNECT LINES AT THE CONNECTION SITE BY TWISTING.**

**5**



**GO TO A SAFE PLACE WITH YOUR EMERGENCY KIT.**

Produced by **BCRenal**  
Provincial Health Services Authority

**fraserhealth**  
Better health. Best in health care.

**Interior Health**

**northern health**  
the northern way of caring

**Providence Health Care**

**Vancouver Coastal Health**  
Promoting wellness. Ensuring care.

**island health**

Provincial Health Services Authority

**BRITISH COLUMBIA**