


Emergency Disconnect Procedure For Patients

for the AK 95S

GRAFT or FISTULA

1



PRESS BLOOD PUMP BUTTON


Blood pump will stop and button will start flashing.

2




USING THE PRODUCTS FOUND IN YOUR CLAMP & CUT KIT...

3



CLOSE TWO CLAMPS ON EACH OF THE TWO BLOODLINES

4



CUT BETWEEN THE CLAMPS

5



GO TO A SAFE PLACE WITH YOUR CLAMP & CUT KIT

Produced by **BCRenal**
an agency of the Provincial Health Services Authority

fraserhealth
Better health. Best in health care.

Interior Health

northern health
the northern way of caring

Providence
HEALTH CARE
How you want to be treated.

Vancouver Coastal Health
Promoting wellness. Ensuring care.

VANCOUVER ISLAND
health authority

Provincial Health Services Authority
Province-wide solutions. Better health.

BRITISH COLUMBIA