


緊急斷離程序

Emergency Disconnect Procedure

型號 AK 98
供病人使用
導管

1




按血泵的按鈕。
血泵將停動，按鈕將開始閃動。

2




使用在急救袋中
找到的用品...

3



關好導管和血路管上的夾子。
在連接位下方和上方加放安全夾。

4



扭開連接位，
將導管和血路管斷離。

5



攜同急救袋
前往安全的地方。

Produced by **BCRenal**
Provincial Health Services Authority

 **fraserhealth**
Better health. Best in health care.

 Interior Health

 **northern health**
the northern way of caring

 **Providence Health Care**

 **Vancouver Coastal Health**
Promoting wellness. Ensuring care.

 **island health**

 Provincial Health Services Authority

 **BRITISH COLUMBIA**